

BoogieBaby™

Nasal Aspirator



BoogieBaby was developed in collaboration with clinicians to include all the options you want in a non-invasive suction device.


- Two sizes: Preemie and Standard
- Flair Tip
- Olive Tip: 2 lengths
- With or without hanging cover

Economical & Easy To Use Oral & Nasal Suction Device & Cover

Manufactured to save money compared to other brands!

The variety of sizes and features make BoogieBaby ideal for use throughout your hospital – NICU, PICU, delivery room, emergency room – wherever suctioning is needed

- Easy-to-use, single-handed suctioning
- Comfortable thumb-control port allows better control for intermittent suctioning
- Preemie and standard tips are soft and flexible, similar to a bulb syringe
- Olive tips minimize the risk of trauma and create a seal at the nares



*Colored covers for visibility.
When used with the convenient hanging cover, BoogieBaby helps comply with VAP protocols.*



How to Order

Item #	Nasal Aspirator Description	Quantity
BB15011100	Preemie tip	Box of 100
BB 15012100	Preemie tip with cover	Box of 100
BB 15013100	Standard tip	Box of 100
BB 15014100	Standard tip with cover	Box of 100
BB 15015100	Flair tip	Box of 100
BB 15016100	Flair tip with cover	Box of 100
BB150112100	Olive tip 4.5 inch	Box of 100
BB150114100	Olive tip 4.5 inch with cover	Box of 100
BB150113100	Olive tip 3.25 inch	Box of 100
BB150115100	Olive tip 3.25 inch with cover	Box of 100

BoogieBaby Preemie & Standard compare to Neotech Little Sucker® Suction Device*

N204	Preemie	Case of 50
N204C	Preemie with cover	Case of 50
N205	Standard	Case of 50
N205C	Standard with cover	Case of 50
N224	Preemie Nasal Tip	Case of 50
N225	Nasal Tip	Case of 50

BoogieBaby Olive Tips compare to Philips BBG® Nasal Aspirator discontinued May 2019**

93032-S	Nasal Aspirator - 4.5" in length	Case of 50
93032-02	Nasal Aspirator - 3.25" in length	Case of 50

Indications for Use

BoogieBaby offers a non-invasive suction alternative for the nasal cavity for use with newborns, infants and children. Indicated for use as a vacuum-powered suction device to aspirate fluids from the mouth and nares.

Instructions For Use

Single Patient Use.

For oral and nasal use only.

Packaged clean & ready to use.

1. Connect to standard suction line.
2. Set suction pressure and use per hospital protocol.
3. Suction sterile water through device after each use.
4. Replace at least every 24 hours.

Product Specifications

Packaged in individual polyethylene bags. Designed & engineered by SandBox Medical in the USA. Manufactured in China at an FDA registered facility.

WARNINGS

Inserting the tip too far into the mouth or nares can cause tissue trauma. The suction device & cover should be kept free of debris. Federal Law (USA) restricts sale of this device to, or on the order of a physician.

*Little Sucker® is a register trademark of Neotech Products. Inc.

**BBG® Nasal Aspirator is a registered trademark of Philips Medical.